



Sensitivity of the sbeadex® forensic kit

L. Forster and E. Johnston, LGC, Queens Road, Teddington, Middlesex, TW11 0LY, UK

Aim

The study set out to evaluate the suitability of the sbeadex® forensic kit for use in extracting DNA from forensic crime stains. The first stage of the evaluation was to test the sensitivity of the sbeadex® forensic kit at recovering DNA from blood and saliva dilution series¹ processed on a liquid handling robotic workstation.

Experimental design

A saliva sample dilution series and blood dilution series were prepared from donors with known genotypes in the concentration range of 2.5 pg/μL to 2000 pg/μL (Table 1).

Each 96-well extraction plate contained 48 saliva samples, 44 blood samples and four negative samples. Two plates were processed using the sbeadex® forensic kit and protocol on a liquid handling robot. The resulting DNA extracts were quantified using the Quantifiler™ Human assay (Applied Biosystems). 25 μL PCRs were performed with 1 ng of DNA added per PCR sample if sufficient levels of DNA were present using SGM Plus® PCR kits (Applied Biosystems). The PCR products were run on an Applied Biosystems 3130xl Genetic Analyser. The STR profiles were then analysed using GeneMapper™ ID v3.2 analysis software.

Saliva sample concentrations	Number of replicates per plate
2.5 pg/μL	8
5 pg/μL	8
10 pg/μL	8
25 pg/μL	8
35 pg/μL	8
50 pg/μL	8
Total	48

Blood sample concentrations	Number of replicates per plate
50 pg/μL	4
100 pg/μL	8
250 pg/μL	8
500 pg/μL	8
750 pg/μL	4
1000 pg/μL	4
1500 pg/μL	4
2000 pg/μL	4
Total	44

Table 1: Details the sample concentrations included in the dilution series

Sample type	Concentration	% of expected alleles observed – plate 1	% of expected alleles observed – plate 2
Saliva	2.5 pg/μL	31	6
	5 pg/μL	53	12
	10 pg/μL	60	25
	25 pg/μL	74	61
	35 pg/μL	82	68
	50 pg/μL	80	63
Blood	50 pg/μL	100	100
	100 pg/μL	100	100
	250 pg/μL	100	100
	500 pg/μL	100	100
	750 pg/μL	100	100
	1000 pg/μL	100	100
	1500 pg/μL	100	100
2000 pg/μL	100	100	

Table 2: Shows the percentage number of expected peaks observed in the STR profiles

Results

Full profiles were obtained for all blood samples across the range of DNA tested 50 pg/μL to 2000 pg/μL. Partial profiles were obtained for all saliva samples across the range of DNA concentration tested 2.5 pg/μL to 50 pg/μL (Table 2).

The mean peak heights for the blood samples varied from 499 rfu to 1726 rfu for the range of DNA tested. The mean peak heights for the saliva samples varied from 39 rfu to 144 rfu for the range of DNA tested (Figure 1).




Figure 1: Shows the average observed peak height by DNA concentration and sample type. Key: s=saliva, b=blood




Figure 3: Example STR profiles generated from saliva samples with DNA concentrations of 5 pg/μL, 10 pg/μL and 50 pg/μL respectively

Conclusion

The investigation has shown that the sbeadex® forensic kit and protocol is able to isolate and purify DNA from samples with a starting concentration of 2.5 pg/μL. Therefore it can be concluded that the sbeadex® forensic kit is suitable for use on forensic crime stain samples where there is only a low level of DNA deposited.




Figure 2: Example STR profiles generated from saliva samples with DNA concentrations of 5 pg/μL, 10 pg/μL and 35 pg/μL respectively




Figure 4: Example STR profiles generated from blood samples with DNA concentrations of 50 pg/μL, 250 pg/μL and 750 pg/μL respectively



LGC Genomics
Germany
Ostendstr. 25 • TGS Haus 8
12459 Berlin

Tel: +49 (0)30 5304 2200
Fax: +49 (0)30 5304 2201
Email: info.de@lgcgenomics.com

United Kingdom
Unit 1-2 Trident Industrial Estate • Pindar Road
Hoddesdon • Herts • EN11 0WZ

Tel: +44 (0) 1992 470757
Fax: +44 (0) 1438 900670
Email: info.uk@lgcgenomics.com

www.lgcgenomics.com

USA
100 Cummings Center • Suite 420H
Beverly • MA 01915

Tel: +1 (978) 232 9430
Fax: +1 (978) 232 9435
Email: info.us@lgcgenomics.com