

ABSTRACT: We introduce a novel probe for qPCR that is embedded with a duplex stabilizing technology and terminated with the Black Hole Quencher®. These BHQ_{plus} probes permit the design of shorter oligonucleotides to detect targets that do not easily accommodate traditional probes. We demonstrate this capability against an AT-rich gene, conserved across several species of closely related bacteria. Fluorescent melting curves reveal the added stabilization of the BHQ_{plus} chemistry and also the consequences of hybridizing to a mismatch, such as the opposite allele of a single-nucleotide polymorphism. Our free PCR application, RealTimeDesign™, is being upgraded to incorporate the thermodynamic parameters specific to this chemistry and support probe-based genotyping. This software program will design numerous SNP assays in a high-throughput fashion, and also accommodate multi-nucleotide polymorphisms. Here, we demonstrate the amplification performance of several SNP assays proposed by RealTimeDesign, with scatter plots resolving the genotypes.

COMPACT

PROBES FOR qPCR: Oligonucleotide sequences for fluorescence-quenched probes are selected to anneal at a precise temperature, usually elevated above that of the primers. BHQ_{plus} probes are synthesized with a duplex-stabilizing chemistry so that truncated oligos continue to anneal at the proper temperature. This feature allows probe design into difficult AT-rich regions.

DUPLEX

STABILIZATION

provided by the BHQ_{plus} chemistry is revealed by ramping the temperature of a mixture containing the probe and its complement, and recording fluorescence as the strands melt apart. Melting curves measured upon mismatched templates reveal single-base discrimination, making BHQ_{plus} probes well-suited for SNP genotyping.

Amplification of an AT-rich target, detected with BHQ_{plus} and traditional probes:





FIGURE 1: Replicate PCR reactions, all amplified using the same primer-set, but detected using either a traditional FAM-BHQ1 probe (black) or a BHQ_{plus} probe (blue). The amplification product is a region within the AT-rich hemolysin III gene of *Bacillus anthracis* (TIGR locus #_BA_2725). The assay is as follows, with the primer-binding regions shown in grey.



At 29 bases in length, the traditional probe remains sufficiently short for FRET quenching, but is predicted to melt from its target below the ideal temperature of 70 °C. In contrast, the BHQ_{plus} probe is only 21 bases in length and yet shows an improved fluorescent response, possibly due to its elevated melting temperature.

Melting curves of FAM- and TET-labeled BHQ_{plus} probes:




FIGURE 2: BHQ_{plus} probes that target a SNP within the gene UGT2B15 were melted in the presence of templates representing each allele. The FAM-labeled probe targets the 'C' allele (blue) while the TET-labeled probe targets the 'T' allele (gold). The melting temperature of each probe, recorded upon their mismatched template is 6-8 °C lower than that of the matched template.

Identical probe sequences were also synthesized without the BHQ_{plus} chemistry, and melted in the presence of their matched templates (black). These melting curves are also shifted toward lower temperatures, confirming the binding stability of BHQ_{plus} probes.

REALTIMEDESIGN

www.qPCRDesign.com

FIGURE 3: Screenshot from RealTimeDesign, displaying the oligo sequences of the genotyping assay targeting UGT2B15, a gene involved in the metabolism of certain drugs such as tamoxifen. Properties of the oligos and amplicons are also presented.

RealTimeDesign is a free software program for qPCR assay design and can be run entirely through a web browser. To assist in the design of genotyping assays, we are developing a module to discriminate SNPs and multi-nucleotide polymorphisms by incorporating the thermodynamics specific to the BHQ_{plus} chemistry. Each allele-specific probe can be labeled with a variety of spectrally-distinct fluorophores so that separate optic channels will record the amplifications of each genotype. Because both probes are present within the PCR reaction, they will compete to bind the amplified template and further enhance the mismatch discrimination.

To confirm the performance of RealTimeDesign and the BHQ_{plus} technology, four distinct SNPs were submitted to its scrutiny, all within genes thought to be involved in drug metabolism. The default assay proposed by the software was then synthesized at Biosearch, with each probe-set labeled with FAM and TET.

FIGURE 4: LEFT: The software avoids certain misalignments between the primers, but presents those that remain. RIGHT: The alignment of the oligos upon the target sequence. Both the sense and antisense strands are considered for optimal placement of the probes.

GENOTYPING

Scatter plots for the four SNP assays record the final fluorescence values from amplification of each allele. All assays were amplified across a seven-point dilution series of the homozygous and heterozygous genotypes, spanning 5,000,000 copies down to 5 copies. The 'FAM' allele is presented in blue, the 'TET' allele is presented in gold, and the heterozygote presented in green.




FIGURE 5: Scatterplot of a SNP assay targeting UGT2B15 allele1 probe: [FAM]-AACATACCTTTCCGTA-[BHQ_{plus}] allele2 probe: [TET]-TTAACATACCTTTTGTA-[BHQ_{plus}]




FIGURE 6: Scatterplot of a SNP assay targeting ABC11 allele1 probe: [FAM]-ATCATCAACCAACGCA-[BHQ_{plus}] allele2 probe: [TET]-ATCATCAACCAACGCG-[BHQ_{plus}]




FIGURE 7: Scatterplot of a SNP assay targeting ATP7B allele1 probe: [FAM]-CCAGACAGAGGTATC-[BHQ_{plus}] allele2 probe: [TET]-CCAGACAGAGGTGTT-[BHQ_{plus}]




FIGURE 8: Scatterplot of a SNP assay targeting CYP2C9 allele1 probe: [FAM]-TGCTTCAGCAGGAGA-[BHQ_{plus}] allele2 probe: [TET]-TGCTTCAGCAGGAGA-[BHQ_{plus}]

METHODS:

All oligo sequences were designed using RealTimeDesign. Synthetic oligonucleotides representing the amplicon sequences were used as templates for each allele. Template dilutions were prepared in water containing yeast tRNA, Roche Molecular systems, with triplicate reactions prepared for all dilution points. All melting curves and PCR amplifications were performed upon a Rotar-Gene 6000, kindly provided by Corbett Research. Melting curve thermal profile: 95 °C for five min, followed by rapid cooling to 40 °C. Temperature was slowly ramped back to 95 °C over the course of three hours with fluorescence recorded every 0.1 °C. PCR thermal cycling profile: 95 °C for five min, followed by 40 cycles of 95 °C for 20 s, 60 °C for 60 s.

Melting curve reaction components:	PCR reaction components:
nuclease-free water	9.00 µL
10X PCR buffer	2.00 µL
MgCl ₂ (50mM)	2.00 µL
Probe (100 µM)	4.00 µL
Template (100 µM)	4.00 µL
Total	20.0 µL
	0.40 µL
	10.0 µL
	5.00 µL
	1.80 µL
	1.80 µL
	0.50 µL
	0.50 µL
Total	20.0 µL

CONCLUSIONS:

The boosted binding stability provided by the BHQ_{plus} chemistry allows these probes to be accommodated into difficult targets and new applications where traditional probes are not well-suited. We demonstrate this capability upon four SNP targets and successfully resolve the genotypes from one another. We have upgraded RealTimeDesign to target single- and multi-nucleotide polymorphisms so that assay design is rapidly automated. BHQ_{plus} may also be suitable for adapting preexisting probe sequences that rely upon other binding stabilizers. We envision other appropriate uses of this technology, including the discrimination of closely-related species and the detection of small targets such as microRNA.

ACKNOWLEDGEMENTS:

We would like to gratefully acknowledge the contributions of all individuals who developed and refined the BHQ_{plus} technology, including Matt Lythe who manufactured the starting material, Darren Dick and Joe Kresch who have optimized the synthesis of this probe type, and Mary Katherine Johanson who provided excellent suggestions and advice. Raymond Peterson, Dean Fiala, and the rest of the team at Celisdon Laboratories are responsible for the sophisticated code that powers RealTimeDesign. info@celisdon-labs.com www.celisdon-labs.com

PCR is a proprietary technology covered by several US patents owned by Roche Molecular Systems, Inc. which have been sub-licensed by PE Corporation in certain fields. Depending on your specific application you may need a license from Roche or PE to practice PCR. Additional information on purchasing licenses to practice the PCR process may be obtained by contacting the Director of Licensing at Roche Molecular Systems, Inc. 1145 Alton Avenue, Alton, CA 94501 or Applied Biosystems business group of the Applied Corporation, 850 Lincoln Center Drive, Foster City, CA 94404. In addition, the 5' nuclease assay and other homogeneous amplification methods used in connection with the PCR process may be covered by U.S. Patents 5,210,015 and 5,487,972, owned by Roche Molecular Systems, Inc. and by U.S. Patent 5,538,840, owned by The Perkin-Elmer Corporation.

Black Hole Quencher, BHQ, and BHQ_{plus} are registered trademarks of Biosearch Technologies, Inc., Novato, CA. Information on licensing programs for the commercial use of these products is available by writing alliances@biosearchtech.com