

-3.372

0.998

100 fa

FIGURE 8. Multiplexed reactions of the CAL Fluor Red 610

false positives or compromise quantitation

cyce FIGURE 10. A small degree of cross-talk is seen bleeding through from CAL Fluor Orange 560 amplifications into the chan nel oriented to detect FAM.

R^2



Independent amplifications of the CAL Fluor Red ay. This pXO1 assay only detects BA Ames and Sterne.





MULTIPLEXING CONSIDERATIONS -

In certain applications such as gene expression measurement, one or more targets may be present in vast excess, creating a challenge for the multiplexed measurement of the other targets. It has been reported that supplementing the master mix with additional reagents, particularly the polymerase, can overcome a loss in sensitivity caused by competing amplifications.



FIGURE 11. LEFT: a diution series of Sterne DNA (1.0 rg to 10 fg) reveals a loss in sensitivity by the CAL Fluor Orange 560 a when there is simultaneous amplification of 1x10⁴⁶ copies of a template representing the pXO2 target (inset). RIGHT: boosting regents in the master mix rescues the amplification of the lower dividuors in competition with this dominant pXO2 amplification [1] K. Persson, K. Hamby, L.A. Ugozzoli, Analytical Biochemistry 344 (2005) 33-42. Jassay boosting the lification CONCLUSIONS

By carefully selecting the target sequences, fluorescent reporters, and the real-time PCR instrument for detection, we have developed a high-performing alTimeDesign was used to rapidly generate TaqMan® assays to the species-specific sequences. To successfully multiplex these assays together, we characterized dye crosstalk, confirmed that each assay has a high amplification efficiency, and identified their detection limit, particularly in the context of disproportionate targets. For situations where one sequence may be in vast excess over the others, we show the benefit of supplementing master mix nents so that depleted reagents don't limit detection sensitivity. Multiplexing demands increased effort to characterize amplification performance but is ideally suited to interrogate multiple genetic signatures from small quantities of sample DNA.

Acknowledgements. We would like to gratefully acknowledge Mary Katherine Johansson, Hans Johansson, and Mark Reddington for their excellent suggestions and advice

FIGURE 8. Multiplexed reactions of the CAL Fluor Red 610 FIGURE 9. Multiplexed reactions of the Quasar 670 assay, am-ssay, amplifying from a ten-fold dilution series of BA Ames DNA plifying from a ten-fold dilution series of BA Ames DNA

- Methods -

3 18

0.997

R^2

Total

Triplicate reactions were prepared for all dilution points and NTCs. Serial dilutions for the construction of standard curves were prepared in nuclease-free water containing 100 ng/µL of yeast tRNA, Roche Molecular systems. Bacillus anthracis Ann and Bacillus anthracis Sterme DNA was obtained through the Multiplexing requires pairing distinct fluorescent reporters that can be resolved from one another. Any crosstalk between channels might cause and Bacillus anthraces Sterme DNA was obtained through the Critical Reagents Program within the Department of Defense. Bacillus coreus strain ATCC 10987 was obtained from ATCC. PCR thermal cycling conditions: S5°C for two min, followed by 60 cycles of: 95°C for 20 s, 60°C for 60 s. All independent reactions were constructed using a master mix with imited reagent concen-trations. All multiplexed reactions were constructed using a master mix with more generous reagent concentuations, except master mix with more generous reagent concentuations, except and an application. Both of these master mix formulations are



20.0 u 20.0 u

"TapMath" is a registered flakfamati of lecone neuroscore "powers" pricelary lechnology on worket by several (50 petiterits owned by Probeh Milecular Systems, inc., which have been sub-licensed by PE Corporation in ordinal fields. Depending on your specificary pointain you may need a lonenet from Robert of PE to patcle CFR Additional Information on pur changing lectroses to patcle the PCR process may be defaued by containing the Director 41. Bowyletins built and a lonenet from Robert of PE to patcle CFR Additional Information on pur changing lectroses to patcle the PCR process may be compared to a long the PCR process Bowyletins built and provide the PCR process star by concentral to Line Centre Dire : Forler Cript CA 94/04, in addition, the 2^e nucleuse says and other homogeneous amplification methods used in mended with the PCR process may be covered by U.S. Pattert 55/08.848, owned by The Perkin-

uencher, BHQ, CAL Fluor, and Quasar are registered trade , Novato, CA. Information on licensing programs for the co vailable by writing alliances@biosearchtech.com